Quick Setup Guide Transmitter M800

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Transmitter M800

Content

1	Safety and Intended Use	37
2	Technical Data and Further Information	41
3	Mounting	41
4	Electrical Connection	41
5	Menu Structure	42
6	Operation On-Site	43
7	Sensor Calibration	45
8	Maintenance	46
9	Disposal	46
10	EC Declaration of Conformity	46

1 Safety and Intended Use

 **Note:** The Quick Setup Guide is a brief operating instruction.

The M800 transmitter must be installed, connected, commissioned, and maintained by qualified specialists e.g. electrician in full compliance with the instructions in this Quick Setup Guide, the applicable norms and legal regulations.

The specialist must have read and understood this Quick Setup Guide and must follow the instructions it contains. If you are unclear on anything in this Quick Setup Guide, you must read the Operation Manual (download from www.mt.com/m800-downloads). The Operation Manual provides detailed information on the instrument.

The M800 transmitter should be operated only by personnel familiar with the transmitter and who are qualified for such work.

Intended Use

The M800 is a 4-wire transmitter for analytical measurements with 4 (0) to 20 mA output signal¹⁾. The M800 is a multi-parameter transmitter and supports the measurements listed in the parameter fit guides. The 2-channel and 4-channel versions are compatible with ISM™ and flow sensors. The 1-channel version is compatible with analog and ISM sensors²⁾.

The M800 transmitter is designed for use in the process industries, in non-hazardous areas.

1) M800 Profinet and Ethernet/IP models do not have analog output

2) M800 Profinet and Ethernet/IP support only ISM sensors, details see Profinet parameter fit guide

M800 parameter fit guide for the 1-channel version

Parameter	Analog	ISM
pH/ORP	•	•
pH/pNa	–	•
UniCond™ 2-e/4-e	– / –	• / •
Conductivity 2-e/4-e	• / •	– / •
Amp. Dissolved Oxygen ppm/ppb/trace	• / • / • ¹⁾	• / • / • ¹⁾
Amp. Oxygen gas ppm/ppb/trace	• / • / • ¹⁾	• / • / • ¹⁾
Optical Dissolved Oxygen	–	• ¹⁾
Dissolved Carbon Dioxide (InPro™ 5000i)	–	•
CO ₂ hi (InPro 5500i)	–	•
Turbidity	• (backscatter)	•

1) Ingold sensors

M800 parameter fit guide for the 2-channel and 4-channel versions

These versions are compatible with the following ISM and flow sensors.

Parameter	Water		Process	
	2-channel	4-channel	2-channel	4-channel
pH/ORP	•	•	•	•
pH/pNa	–	–	•	•
UniCond 2-e	•	•	•	•
Conductivity 4-e	•	•	•	•
Amp. Dissolved Oxygen ppm/ppb/trace	– / • / – ²⁾	– / • / – ²⁾	• / • / • ¹⁾	• / • / • ¹⁾
Amp. Oxygen gas ppm/ppb/trace	– / • / – ²⁾	– / • / – ²⁾	• / • / • ¹⁾	• / • / • ¹⁾
Optical Dissolved Oxygen	• ²⁾	• ²⁾	• ^{1), 3)}	• ^{1), 3)}
Dissolved Carbon Dioxide (InPro 5000i)	–	–	•	•
CO ₂ hi (InPro 5500i)	–	–	• ³⁾	• ³⁾
TOC	•	•	–	–
Dissolved Ozone	•	•	–	–
Flow	•	•	–	–

1) Ingold sensors

2) Thorlon sensors

3) 2-channel: An opt. dissolved sensor or a CO₂ hi sensor has to be connected to channel 2. 4-channel: Optical dissolved sensors and CO₂ hi sensors have to be connected to channel 2 and/or to channel 4

M800 Profinet and Ethernet/IP parameter fit guide

Parameter	Water		Process	
	2-channel	4-channel	1-channel	2-channel
pH/ORP	•	•	•	•
pH/pNa	–	–	•	•
UniCond™ 2-e/4-e	•	•	•	•
Conductivity 4-e	•	•	•	•
Amp. Dissolved Oxygen ppm/ppb/trace	– / • / – ²⁾	– / • / – ²⁾	• / • / • ¹⁾	• / • / • ¹⁾
Amp. Oxygen gas ppm/ppb/trace	– / • / – ²⁾	– / • / – ²⁾	• / • / • ¹⁾	• / • / • ¹⁾
Optical Dissolved Oxygen	• ^{2),3)}	• ²⁾	• ^{1),3)}	• ^{1),3)}
Dissolved Carbon Dioxide (InPro 5000i)	–	–	•	•
CO ₂ hi (InPro 5500i)	–	–	• ³⁾	• ³⁾
InPro 86X0i	–	–	•	•
TOC	•	•	–	–
Dissolved Ozone	•	•	–	–
Flow	•	•	–	–

1) Ingold sensors

2) Thornton sensors

3) Only one optical DO or thermal conductivity CO₂ sensor can be used together with 2-channel transmitter

M800 21CFR RecordLOC™ parameter fit guide

Parameter	Water
	2-channel
pH/ORP	–
pH/pNa	–
UniCond 2-e/4-e	•
Conductivity 4-e	–
Amp. Dissolved Oxygen ppm/ppb/trace	–
Amp. Oxygen gas ppm/ppb/trace	–
Optical Dissolved Oxygen	–
Dissolved Carbon Dioxide (InPro 5000i)	–
CO ₂ hi (InPro 5500i)	–
InPro 86X0i	–
TOC	•
Dissolved Ozone	• ¹⁾
Flow	–

1) Only ISM sensors

2 Technical Data and Further Information

The most relevant technical data such as supply voltage are shown on the name plate either outside or inside the transmitter housing. For further technical data such as accuracy refer to the Operation Manual. This document, the Operation Manual and software can be downloaded from www.mf.com/m800-downloads.

3 Mounting

The M800 transmitter is available as ½ DIN version.
For installation drawings refer to the Operation Manual.

DANGER! Mortal danger by electric shock or risk of electrical shock: The maximum screw-in depth of the mounting holes in the housing is 12 mm (0.47 inch). Do not exceed maximum screw-in depth.

1. Mount the supplied cable glands at the housing.
2. Mount the transmitter. You have the following possibilities:
 - Panel mounting, wall mounting or pipe mounting

4 Electrical Connection


DANGER! Mortal danger by electric shock: Power off instrument during electrical connection.

Note: Do not supply power to the **Aout** terminals.

1. Switch off supply voltage.
2. Connect mains supply as follows:
 - 20 to 30 V DC: **N (-)** for Neutral and **L (+)** for Line
 - 100 to 240 V AC: **N** for Neutral and **L** for Line
3. Connect sensor, analog output signals, digital input signals and relay output signals according to the Operation Manual.

5 Menu Structure


1) For 21CFR RecordLOC only

2) For Ethernet/IP Setting

3) For Profinet only

6 Operation On-Site

A Start screen (example)

- 1 Changing between the channels, only 2-channel and 4-channel versions
- 2 1st line, standard configuration
- 3 2nd line, standard configuration
- 4 3rd line, depends on configuration
- 5 4th line, depends on configuration
- 6 Profinet or Ethernet/IP connection indicator¹⁾

B Menu screen (example)

C ISM Menu screen

- 1) For Profinet and Ethernet/IP models

Note: Above screen(example) may vary based on specific models.
Pls refer to required model's installation manual for detailed information.

Operating element	Description
	Enter Menu screen
	Enter Start screen
ISM	Enter ISM menu
	Enter Favorite menu
	Enter Calibration menu
	Enter Configuration menu
	Return to Menu screen
	Enter next-lower menu level, here e.g. iMonitor, Messages or ISM Diagnostics
	Return to next-higher menu level
	<ul style="list-style-type: none"> • Change between pages within one menu level • 2-channel and 4-channel: Change between the channels

7 Sensor Calibration

 **Note:** For best process calibration results observe the following points. Take sample as close as possible to the measurement point of the sensor. Measure the sample at process temperature.

For a detailed description of the “Process”, “1-Point” and “2-Point” calibration method see Operation Manual of the M800 transmitter. Via iSense software you can calibrate the sensor with the “1-Point” or “2-Point” calibrate method. See Operation Manual of the iSense software.

As soon as the calibration is in progress no other calibration can be started.

Sensor Calibration Menu

After every successful calibration different options are available. If “Adjust”, “SaveCal” or “Calibrate” is chosen, the message “Calibration saved successfully! Reinstall sensor” is displayed. Press “Done”.

Option	Analog sensors (1-channel version only)	ISM (digital) sensors
Analog sensors: Save Cal	Calibration values are stored in the transmitter and used for the measurement. Additionally, the calibration values are stored in the calibration data.	Calibration values are stored in the sensor and used for computing measurements. Additionally, the calibration values are stored in the calibration history.
ISM sensors: Adjust		
Calibrate	The function “Calibrate” is not applicable for analog sensors.	Calibration values are stored in the calibration history for documentation are not used for computing measurements. The calibration values from the last adjustment are still used for computing measurements.
Cancel	Calibration values are discarded.	Calibration values are discarded.

8 Maintenance

The transmitter requires no maintenance.

Clean the surfaces with a soft damp cloth and dry the surfaces with a cloth carefully.

9 Disposal

Observe the applicable local or national regulations concerning the disposal of "Waste electrical and electronic equipment".

Disassemble the transmitter according to resources. Sort resources and supply them to recycling. Non-recyclable materials have to be disposed of in an environmental-friendly manner.

10 EC Declaration of Conformity

The EC Declaration of conformity is part of the delivery.

ISM, InPro and UniCond are trademarks of the METTLER TOLEDO Group.
All other trademarks are the property of their respective owners.

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© 06/2023 METTLER TOLEDO. All rights reserved
Mettler-Toledo GmbH
Mettler-Toledo Thornton, Inc.